

**产品名称:** Postmeiotic Segregation Increased 2(PMS2) Mouse Monoclonal Antibody  
**产品货号:** AMM22040

### 产品概述 (Summary)

|                   |                                                                                                      |
|-------------------|------------------------------------------------------------------------------------------------------|
| 描述(Description)   | Mouse Monoclonal Antibody                                                                            |
| 宿主(Host)          | Mouse                                                                                                |
| 应用(Application)   | IHC,ELISA                                                                                            |
| 种属反应(Reactivity)  | Human,Mouse                                                                                          |
| 偶联物(Conjugation)  | Unconjugated                                                                                         |
| 修饰(Modification)  | Unmodified                                                                                           |
| 同种型(Isotype)      | IgG1,Kappa                                                                                           |
| 克隆(Clonality)     | Monoclonal                                                                                           |
| 剂型(Form)          | Liquid                                                                                               |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液(Buffer)      | PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA                                                       |
| 纯化 (Purification) | The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen. |

### 产品应用(Application)

|                       |                                        |
|-----------------------|----------------------------------------|
| 稀释比(Dilution Ratio)   | IHC 1:50-100;ELISA 1:500-5000          |
| 分子量(Molecular Weight) | Calculated MW:96kDa,Observed MW:110kDa |

### 抗原信息(Antigen Information)

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| 基因名(Gene Name)        | PMS2 PMSL2                                                                                         |
| 别名(Alternative Names) | Mismatch repair endonuclease PMS2;DNA mismatch repair protein PMS2;PMS1 protein homolog 2;         |
| 基因 ID(Gene ID)        | Human:5395                                                                                         |
| SwissProt ID          | Human:P54278                                                                                       |
| 免疫原(Immunogen)        | Synthesized peptide derived from human Postmeiotic Segregation Increased 2(PMS2) AA range: 600-700 |

### 研究背景 (Background)

The protein encoded by this gene is a key component of the mismatch repair system that functions to correct DNA mismatches

and small insertions and deletions that can occur during DNA replication and homologous recombination. This protein forms heterodimers with the gene product of the mutL homolog 1 (MLH1) gene to form the MutL-alpha heterodimer. The MutL-alpha heterodimer possesses an endonucleolytic activity that is activated following recognition of mismatches and insertion/deletion loops by the MutS-alpha and MutS-beta heterodimers, and is necessary for removal of the mismatched DNA. There is a DQHA(X)2E(X)4E motif found at the C-terminus of the protein encoded by this gene that forms part of the active site of the nuclease. Mutations in this gene have been associated with hereditary nonpolyposis colorectal cancer (HNPCC; also known as Lynch syndrome) and Turcot syndrome


### 研究领域 (Research Area)

Pathology

### 注意事项 (Note)

For research use only.

### 图片 (Image Data)



Human colon carcinoma tissue was stained with anti-PMS2 Antibody

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>