

产品名称: Ki-67 Rabbit Monoclonal Antibody

产品货号: AMRe36647

产品概述 (Summary)

描述(Description)	Rabbit Monoclonal Antibody
宿主(Host)	Rabbit
应用(Application)	IHC, WB, ELISA
种属反应(Reactivity)	Human
偶联物(Conjugation)	Unconjugated
修饰(Modification)	Unmodified
同种型(Isotype)	IgG1, Kappa
克隆(Clonality)	Monoclonal
剂型(Form)	Liquid
保存条件(Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液(Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% BSA
纯化 (Purification)	Recombinant Expression and Affinity purified

产品应用(Application)

稀释比(Dilution Ratio)	IHC 1:100-500; WB 1:500-1000; ELISA 1:5000-20000
分子量(Molecular Weight)	Calculated MW: 359kDa

抗原信息(Antigen Information)

基因名(Gene Name)	MKI67
别名(Alternative Names)	Antigen identified by monoclonal Ki 67; Antigen identified by monoclonal Ki-67; Antigen KI-67; Antigen KI67; Antigen KI67; KI67_HUMAN; KIA; Marker of proliferation Ki-67; MIB1; MIB; MKI67; PPP1R105; Proliferation marker protein Ki-67; Proliferation related Ki 67 antigen; Protein phosphatase 1 regulatory subunit 105; RP11-380J17.2
基因 ID(Gene ID)	Human: 4288
SwissProt ID	Human: P46013
免疫原(Immunogen)	Synthesized peptide derived from human Ki-67 AA range: 1200-1400

研究背景 (Background)

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced

transcript variants have been described. A related pseudogene exists on chromosome X. [provided by RefSeq, Mar 2009],

研究领域 (Research Area)

Pathology

注意事项 (Note)

For research use only.

图片 (Image Data)

Human tonsil tissue was stained with Anti-Ki-67 Antibody

EnkiLife 优势产品推荐

产品类别	产品货号	产品名称
WB 解决方案	RA10020	2h 极速 WB 即用型全流程试剂盒
	RA10021	4h 快速 WB 即用型全流程试剂盒
	RA10042	5h 畅享版 WB 全流程试剂盒
	RA10037	校准级彩色预染蛋白 Marker (8-180kDa)
	RA10038	校准级彩色预染蛋白 Marker (10-250kDa)
	RA10039	校准级高分子彩色预染蛋白 Marker (25-400kDa)
TSA 多重荧光染色试剂盒	RA10008	TSA 双标三色多重荧光染色试剂盒 (mIHC)
	RA10009	TSA 三标四色多重荧光染色试剂盒 (mIHC)
	RA10010	TSA 四标五色多重荧光染色试剂盒 (mIHC)
	RA10011	TSA 五标六色多重荧光染色试剂盒 (mIHC)
	RA10012	TSA 六标七色多重荧光染色试剂盒 (mIHC)
IHC 检测试剂盒	RA10006	HRP Anti-Mouse/Rabbit IHC Detection System
	RA10007	Polymer-HRP Anti-Mouse/Rabbit IHC Detection System
抗体标记试剂盒	RE80004p	辣根过氧化物酶(HRP)抗体标记试剂盒
	RE80002q	Sulfo-NHS-生物素标记试剂盒
	RE80007p	Cy3 荧光素标记试剂盒
	RE80011p	Fluor488 荧光素标记试剂盒
	RE80017p	Fluor750 荧光素标记试剂盒
	RE80005p	藻红蛋白(R-PE) 抗体快速标记试剂盒
	RE80040	PE-Cy7 串联染料抗体快速标记试剂盒
稳转细胞系构建服务 (免费赠送全膜 WB 验证)	TS-0001	过表达稳转细胞系构建
	TS-0002	敲低稳转细胞系构建
	TS-0003	敲除细胞系构建

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>