

产品名称: FGFR4 Mouse Monoclonal Antibody
产品货号: AMM00718



产品概述 (Summary)

产品名称 (Production Name)	FGFR4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Ascitic Fluid

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGFR4
别名 (Alternative Names)	FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD334
基因 ID (Gene ID)	2264
蛋白 ID (SwissProt ID)	P22455.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Receptor for acidic fibroblast growth factor. Does not bind to basic fibroblast growth factor. Binds FGF19.


产品名称: FGFR4 Mouse Monoclonal Antibody
产品货号: AMM00718




研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of FGFR4 in Human FGFR4 (aa22369) lysates using FGFR4 antibody.



Immunofluorescence analysis of FGFR4 in HEK293 cells transfected with FGFR4hlgGfc using FGFR4 antibody (green) Blue: DRAQ5 fluorescent DNA dye.

注意事项 (Note)

For research use only .