

产品名称: Ki67 (10B6) Mouse Monoclonal Antibody  
产品货号: AMM00719



## 产品概述 (Summary)

|                        |                                       |
|------------------------|---------------------------------------|
| 产品名称 (Production Name) | Ki67 (10B6) Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody             |
| 宿主 (Host)              | Mouse                                 |
| 应用 (Application)       | IHC                                   |
| 种属反应性 (Reactivity)     | Human                                 |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                     |
| 同种型 (Isotype)       | IgG2b                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Ascitic Fluid                                                                                  |

## 免疫原信息 (Immunogen)

|                        |                      |
|------------------------|----------------------|
| 基因名 (Gene Name)        | MKI67                |
| 别名 (Alternative Names) | MKI67; Antigen KI-67 |
| 基因 ID (Gene ID)        | 4288                 |
| 蛋白 ID (SwissProt ID)   | P46013.              |

## 产品应用 (Application)

|                          |                |
|--------------------------|----------------|
| 稀释比 (Dilution Ratio)     | IHC 1:50-1:100 |
| 蛋白分子量 (Molecular Weight) | -              |

## 研究背景 (Background)

Required to maintain individual mitotic chromosomes dispersed in the cytoplasm following nuclear envelope

**产品名称: Ki67 (10B6) Mouse Monoclonal Antibody**  
**产品货号: AMM00719**




disassembly(PubMed:27362226). Associates with the surface of the mitotic chromosome, the perichromosomal layer, and covers a substantial fraction of the chromosome surface (PubMed:27362226).

## 研究领域 (Research Area)

Cell Biology

## 图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human lymph node (A), esophagus (B), lung cancer (C), rectum cancer (D), using Ki67 (10B6) antibody with DAB staining. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## 注意事项 (Note)

For research use only .