

产品名称: NCK1 (7H2) Mouse Monoclonal Antibody
产品货号: AMM00720



产品概述 (Summary)

产品名称 (Production Name)	NCK1 (7H2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,FC
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCK1
别名 (Alternative Names)	NCK; nck-1; NCKalpha
基因 ID (Gene ID)	4690
蛋白 ID (SwissProt ID)	P16333.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Adapter protein which associates with tyrosine-phosphorylated growth factor receptors or their cellular substrates.

产品名称: NCK1 (7H2) Mouse Monoclonal Antibody
产品货号: AMM00720




Maintains low levels of EIF2S1 phosphorylation by promoting its dephosphorylation by PP1.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of NCK1 (7H2) in Human NCK1 (AA: 20337) recombinant protein.) using NCK1 (7H2) antibody.




Western blot analysis of NCK1 (7H2) in Jurkat , HeLa , HEK293 , A431 , K562 ,and COS7 lysates using NCK1 (7H2) antibody



Immunohistochemistry analysis of paraffin-embedded rectum cancer tissues using NCK1 antibody with DAB staining.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: NCK1 (7H2) Mouse Monoclonal Antibody
产品货号: AMM00720



Flow cytometry analysis of Jurkat stained with NCK1 antibody (green) and negative control (red).

注意事项 (Note)

For research use only .