

产品名称: alpha Actinin (1A5) Mouse Monoclonal Antibody
产品货号: AMM00737



产品概述 (Summary)

产品名称 (Production Name)	alpha Actinin (1A5) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTN1
别名 (Alternative Names)	actinin 1 smooth muscle; actinin alpha 1; ACTN 1; FLJ40884
基因 ID (Gene ID)	87
蛋白 ID (SwissProt ID)	P12814.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling

产品名称: alpha Actinin (1A5) Mouse Monoclonal Antibody
产品货号: AMM00737




protein.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using alpha Actinin (1A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using alpha Actinin (1A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .