

产品名称: Cyclin B1 (1B3) Mouse Monoclonal Antibody
产品货号: AMM00754



产品概述 (Summary)

产品名称 (Production Name)	Cyclin B1 (1B3) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCNB1
别名 (Alternative Names)	CCNB1; CCNB; G2/mitotic-specific cyclin-B1
基因 ID (Gene ID)	891
蛋白 ID (SwissProt ID)	P14635.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Essential for the control of the cell cycle at the G2/M (mitosis) transition. a member of the highly conserved cyclin family,

产品名称: Cyclin B1 (1B3) Mouse Monoclonal Antibody
产品货号: AMM00754




whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using Cyclin B1 (1B3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Cyclin B1 (1B3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .