

产品名称: Dynamin 1 (5B9) Mouse Monoclonal Antibody
产品货号: AMM00857



产品概述 (Summary)

产品名称 (Production Name)	Dynamin 1 (5B9) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Ascitic Fluid

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNM1
别名 (Alternative Names)	DNM1; DNM; Dynamin-1
基因 ID (Gene ID)	1759
蛋白 ID (SwissProt ID)	Q05193.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 97 kDa; Observed MW: 97 kDa

研究背景 (Background)

Microtubule-associated force-producing protein involved in producing microtubule bundles and able to bind and

产品名称: **Dynamin 1 (5B9) Mouse Monoclonal Antibody**
产品货号: **AMM00857**




hydrolyze GTP. Most probably involved in vesicular trafficking processes. Involved in receptor-mediated endocytosis.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Dynamin 1 (5B9) in C6 , NIH/3T3 , SK-N-SH , LN18 , SH-SY5Y lysates and rat brain lysates using Dynamin 1 (5B9) antibody.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Dynamin 1 (5B9) antibody, with DAB staining. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .