

产品名称: ApoA4 Mouse Monoclonal Antibody
产品货号: AMM00861



产品概述 (Summary)

产品名称 (Production Name)	ApoA4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Ascitic Fluid

免疫原信息 (Immunogen)

基因名 (Gene Name)	APOA4
别名 (Alternative Names)	Apo AIV; APOA 4; ApoA IV; APOA4; Apolipoprotein A IV; Apolipoprotein A4
基因 ID (Gene ID)	337
蛋白 ID (SwissProt ID)	P06727.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)

APOA4 (apolipoprotein A-IV) is a component of HDL and chylomicrons. Its primary site of synthesis is the intestine, in

产品名称: ApoA4 Mouse Monoclonal Antibody
产品货号: AMM00861




association with lymph chylomicron particles.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of ApoA4 in Human serum , Human plasma , HepG2 lysates and SMMC7721 lysates using ApoA4 antibody.

注意事项 (Note)

For research use only .