

产品名称: Cytokeratin 6 (1A5) Mouse Monoclonal Antibody
产品货号: AMM00865



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 6 (1A5) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT6A/KRT6B/KRT6C CK 6A; CK 6B; CK 6C; CK 6D; CK 6E; CK-6C; CK-6E; K6a keratin; K6b keratin; K6C; K6c keratin; K6d keratin; K6e keratin; Keratin; KRT6A; KRT6B; KRT6C; KRT6D; KRT6E
别名 (Alternative Names)	
基因 ID (Gene ID)	3853/3854/286887
蛋白 ID (SwissProt ID)	P02538/P04259/P48668.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa; Observed MW: 56-60 kDa

产品名称: Cytokeratin 6 (1A5) Mouse Monoclonal Antibody
产品货号: AMM00865


研究背景 (Background)

There are at least six isoforms of human type II keratin-6 (K6). There are two types of cytoskeletal and microfibrillar keratin, I (acidic) and II (neutral to basic) (40-55 and 56-70 kDa, respectively).


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded Human breastcancer tissue using Cytokeratin 6 (1A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of Cytokeratin 6 (1A5) in Human breastcancer tissue using Cytokeratin 6 (1A5) antibody (red), and DAPI (blue).

产品名称: Cytokeratin 6 (1A5) Mouse Monoclonal Antibody
产品货号: AMM00865



Western blot analysis of Cytokeratin 6 (1A5) in Hela, C6, 3T3, K562, HepG2 lysates using Cytokeratin 6 (1A5) antibody

注意事项 (Note)

For research use only .