

产品名称: Bovine Serum Albumin (10G8) Mouse Monoclonal Antibody
产品货号: AMM00886



产品概述 (Summary)

产品名称 (Production Name)	Bovine Serum Albumin (10G8) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Bovine

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Alb
别名 (Alternative Names)	BSA; ALB; Albumin; HSA; PRO0883; Serum Albumin
基因 ID (Gene ID)	280717.0
蛋白 ID (SwissProt ID)	P02769.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa; Observed MW: 69 kDa

研究背景 (Background)

产品名称: Bovine Serum Albumin (10G8) Mouse Monoclonal Antibody
产品货号: AMM00886




Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca²⁺, Na⁺, K⁺, fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Bovine Serum Albumin (10G8) in Bovine Serum Albuminlysates using Bovine Serum Albumin (10G8) antibody.

注意事项 (Note)

For research use only .