

产品名称: Ferritin Light Chain (1A2) Mouse Monoclonal Antibody
产品货号: AMM03303



产品概述 (Summary)

产品名称 (Production Name)	Ferritin Light Chain (1A2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FTL
别名 (Alternative Names)	Ferritin L chain; Ferritin L subunit; Ferritin light chain; Ferritin light polypeptide; Ferritin light polypeptide like 3; FRIL; FRIL_HUMAN; FTL; L apoferritin; LFTD; MGC71996; NBIA 3; NBIA3.
基因 ID (Gene ID)	2512
蛋白 ID (SwissProt ID)	P02792.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

产品名称: Ferritin Light Chain (1A2) Mouse Monoclonal Antibody
产品货号: AMM03303




研究背景 (Background)

Stores iron in a soluble, non-toxic, readily available form. Important for iron homeostasis. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. Also plays a role in delivery of iron to cells. Mediates iron uptake in capsule cells of the developing kidney.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Indirect ELISA assay for antiFerritin Light Chain mouse mAb. Antigen coating concentration: 4ug/ml.

注意事项 (Note)

For research use only .