

产品名称: Anthrax Protective gen (5C1) Mouse Monoclonal Antibody
产品货号: AMM03305



产品概述 (Summary)

产品名称 (Production Name)	Anthrax Protective gen (5C1) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Bacillus anthracis

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	pagA
别名 (Alternative Names)	Protective antigen
基因 ID (Gene ID)	45025512.0
蛋白 ID (SwissProt ID)	P13423.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

产品名称: Anthrax Protective gen (5C1) Mouse Monoclonal Antibody
产品货号: AMM03305




Protective antigen constitutes one of the three proteins composing the anthrax toxin; it mediates attachment to host cells and translocation of edema factor (EF) and lethal factor (LF) into the host cytoplasm.

研究领域 (Research Area)

Microbiology

图片 (Image Data)



Standard Curve for Anthrax Protective Antigen: Capture Antibody mouse mAb 201240 to Anthrax Protective Antigen at 4ug/ml and other antibody for detecting.

注意事项 (Note)

For research use only .