

产品名称: FABP3 (3E8B10) Mouse Monoclonal Antibody
产品货号: AMM03312



产品概述 (Summary)

产品名称 (Production Name)	FABP3 (3E8B10) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FABP3 422 protein; Cardiac Fatty Acid Binding Protein; FABP 11; FABP 3; FABP11; FABP3; FABPH_HUMAN; fatty acid binding protein 11; Fatty acid binding protein 3; Fatty acid binding protein 3 muscle and heart; Fatty acid binding protein 3 muscle and heart mammary derived growth inhibitor; Fatty acid binding protein 3 muscle; Fatty acid binding protein 3
别名 (Alternative Names)	
基因 ID (Gene ID)	2170
蛋白 ID (SwissProt ID)	P05413.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
----------------------	----------------------

产品名称: FABP3 (3E8B10) Mouse Monoclonal Antibody
产品货号: AMM03312



蛋白分子量 (Molecular Weight) -


研究背景 (Background)

FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Standard Curve for HFABP: Capture Antibody mouse mAb (5C4F5) to HFABP at 4ug/ml and Detector Antibody mouse mAb(3E8B10u30017B10G11u30014H1B8)to HFABP at 0.08 ug/ml.

注意事项 (Note)

For research use only .