

产品名称: beta 2 Microglobulin (4C3) Mouse Monoclonal Antibody
产品货号: AMM03314



产品概述 (Summary)

产品名称 (Production Name)	beta 2 Microglobulin (4C3) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	B2M
别名 (Alternative Names)	B2MG; Beta 2 microglobulin; Beta 2 microglobulin; Beta-2-microglobulin form pl 5.3; CDABP0092; Hdcma22p
基因 ID (Gene ID)	567
蛋白 ID (SwissProt ID)	P61769.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

产品名称: beta 2 Microglobulin (4C3) Mouse Monoclonal Antibody
产品货号: AMM03314




研究背景 (Background)

Major histocompatibility complex (MHC) class 1 molecules bind to antigens for presentation on the surface of cells. The proteasome is responsible for producing these antigens from the components of foreign pathogens. MHC class 1 molecules consist of an a heavy chain that contains three subdomains ($\alpha 1$, $\alpha 2$, $\alpha 3$), and a non-covalent associating light chain, known as β -2-Microglobulin.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



sELISA standard Curve for u03b22microglobulin.

注意事项 (Note)

For research use only .