

产品名称: RBP4 (8C12) Mouse Monoclonal Antibody
产品货号: AMM03323



产品概述 (Summary)

产品名称 (Production Name)	RBP4 (8C12) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBP4 OTTHUMP00000020114; OTTHUMP00000020115; OTTHUMP00000020116; Plasma retinol binding protein 4; Plasma retinol-binding protein; Plasma retinol-binding protein(1-176); PRBP; PRO2222; RBP; RBP4; RDCCAS; RET4_HUMAN; Retinol binding protein 4; Retinol binding protein 4 interstitial; Retinol binding protein 4 plasma.
别名 (Alternative Names)	
基因 ID (Gene ID)	5950
蛋白 ID (SwissProt ID)	P02753.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
----------------------	----------------------

产品名称: RBP4 (8C12) Mouse Monoclonal Antibody
产品货号: AMM03323



蛋白分子量 (Molecular Weight) -


研究背景 (Background)

Delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin, this prevents its loss by filtration through the kidney glomeruli.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Standard Curve for RBP4: Capture Antibody mouse mAb (8C12C11) to RBP4 and Detector Antibody mouse mAb(7C5F6) to RBP4.

注意事项 (Note)

For research use only .