

产品名称: FABP1 (8G6) Mouse Monoclonal Antibody
产品货号: AMM03328



产品概述 (Summary)

产品名称 (Production Name)	FABP1 (8G6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FABP1 FABP 1; FABP1; FABP1; FABPL; FABPL_HUMAN; Fatty Acid Binding Protein 1; Fatty acid binding protein 1 liver; Fatty Acid Binding Protein; Fatty acid-binding protein 1; Fatty acid-binding protein; Fatty acid-binding protein liver; L FABP; L-FABP; liver; Liver-type fatty acid-binding protein.
别名 (Alternative Names)	
基因 ID (Gene ID)	2168
蛋白 ID (SwissProt ID)	P07148.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

产品名称: FABP1 (8G6) Mouse Monoclonal Antibody
产品货号: AMM03328




研究背景 (Background)

Binds free fatty acids and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Standard Curve for LFABP: Capture Antibody mouse mAb (8G6D5) to LFABP and Detector Antibody mouse mAb(2D6F6) to LFABP.



Observed LFABP levels in several Human urine samples and Human serum samples at different dilution.

注意事项 (Note)

For research use only .