

产品名称: C Reactive Protein (1H3) Mouse Monoclonal Antibody
产品货号: AMM03329



产品概述 (Summary)

产品名称 (Production Name)	C Reactive Protein (1H3) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRP
别名 (Alternative Names)	CRP; PTX1; C-reactive protein
基因 ID (Gene ID)	1401
蛋白 ID (SwissProt ID)	P02741.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

产品名称: C Reactive Protein (1H3) Mouse Monoclonal Antibody
产品货号: AMM03329




Serum CRP values are routinely measured to detect or monitor many human diseases and minor elevations of C-reactive protein are predictive of cardiovascular events in patients with coronary heart disease. CRP promotes host defense, is an efficient activator of the classical C cascade, and shown to be opsonic. The opsonic and C activation activities are responsible for the nonspecific host defense activities providing protection in the preimmune stage of acute infections.


研究领域 (Research Area)

Immunology

图片 (Image Data)



sELISA standard Curve for CReactive Protein: Capture Antibody mouse mAb to CReactive Protein at 1ug/ml and 2012781H3 was used for detecting.



sELISA standard Curve for CReactive Protein: Capture Antibody Rabbit pAb to CReactive Protein at 4ug/ml and 2012781H3 was used for detecting.

注意事项 (Note)

For research use only .