

产品名称: Natriuretic Peptides B (3C9) Mouse Monoclonal Antibody
产品货号: AMM03332



产品概述 (Summary)

产品名称 (Production Name)	Natriuretic Peptides B (3C9) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPPB
别名 (Alternative Names)	NPPB; Natriuretic peptides B; Gamma-brain natriuretic peptide
基因 ID (Gene ID)	4879
蛋白 ID (SwissProt ID)	P16860.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)


产品名称: Natriuretic Peptides B (3C9) Mouse Monoclonal Antibody
产品货号: AMM03332

Brain natriuretic peptide (BNP) circulates in blood as a peptide hormone with natriuretic, vasodilatory and renin inhibitory properties. BNP is secreted predominantly by the left ventricular myocytes in response to volume expansion and pressure overload.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



sELISA standard Curve for NTproBNP: Capture Antibody mouse mAb (2012893C9) to NTproBNP and Detector Antibody mouse mAb(2012891B5) to NTproBNP.



AntiFITC mouse mAb was used to capture FITClabeled NTproBNP mouse mAb (2012893C9), HRPlabeled Natriuretic Peptides B mouse mAb (2012891B5) was used for detecting.

注意事项 (Note)

For research use only .