

产品名称: Human IgM (3E12) Mouse Monoclonal Antibody
产品货号: AMM03336



产品概述 (Summary)

产品名称 (Production Name)	Human IgM (3E12) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGHM
别名 (Alternative Names)	Ig mu chain C region
基因 ID (Gene ID)	3507.0
蛋白 ID (SwissProt ID)	P01871.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

IgM measurement yields information about the body's immediate resistance and response to infection as well as

产品名称: Human IgM (3E12) Mouse Monoclonal Antibody
产品货号: AMM03336




information related to specific diseases. Decreased levels are associated with immune deficiency states, hereditary deficiencies, and myeloma. Increased levels can be associated with Waldenstrom's macroglobulinemia, chronic infection and hepatocellular disease.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Indirect ELISA assay for mouse Human IgM antibody. Antigen coating concentration: 2ug/ml.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Human IgM (3E12) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .