

产品名称: GFAP (9A2) Mouse Monoclonal Antibody
产品货号: AMM03339



产品概述 (Summary)

产品名称 (Production Name)	GFAP (9A2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GFAP
别名 (Alternative Names)	GFAP; FLJ45472; cb345; ALXDRD
基因 ID (Gene ID)	2670
蛋白 ID (SwissProt ID)	P14136.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

研究背景 (Background)

GFAP is commonly used as a marker for intracranial and intraspinal tumors arising from astrocytes. In addition, GFAP

产品名称: GFAP (9A2) Mouse Monoclonal Antibody
产品货号: AMM03339




intermediate filaments are also present in nonmyelin-forming Schwann cells in the peripheral nervous system


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human liver tissue using GFAP (9A2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of GFAP (9A2) in mouse brain tissue using GFAP antibody(5C8)(red), and DAPI (blue).

注意事项 (Note)

For research use only .