

产品名称: Phospho-Tyrosine (6G6) Mouse Monoclonal Antibody
产品货号: AMM03342



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Tyrosine (6G6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	-
别名 (Alternative Names)	Phospho-Tyrosine
基因 ID (Gene ID)	-
蛋白 ID (SwissProt ID)	-.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Refer to figures

研究背景 (Background)

产品名称: Phospho-Tyrosine (6G6) Mouse Monoclonal Antibody
产品货号: AMM03342




Tyrosine phosphorylation plays a key role in cellular signaling. Research studies have shown that in cancer, unregulated tyrosine kinase activity can drive malignancy and tumor formation by generating inappropriate proliferation and survival signals. Antibodies specific for phospho-tyrosine have been invaluable reagents in these studies.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Phospho-Tyrosine (6G6) in Jurkat, 293T, rat Brain, mouse Brain lysates using Phospho-tyrosine antibody.



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma using Phospho-tyrosine antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .