

产品名称: TBL1X (4H2) Mouse Monoclonal Antibody
产品货号: AMM03414



产品概述 (Summary)

产品名称 (Production Name)	TBL1X (4H2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBL1X
别名 (Alternative Names)	TBL1X; TBL1; F-box-like/WD repeat-containing protein TBL1X; SMAP55; Transducin beta-like protein 1X; Transducin-beta-like protein 1; X-linked
基因 ID (Gene ID)	6907
蛋白 ID (SwissProt ID)	O60907.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 58 kDa

研究背景 (Background)

产品名称: TBL1X (4H2) Mouse Monoclonal Antibody
产品货号: AMM03414




F-box-like protein involved in the recruitment of the ubiquitin/19S proteasome complex to nuclear receptor-regulated transcription units.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of TBL1X in SW480 lysates using TBL1X antibody.



Immunocytochemistry analysis of TBL1X (4H2) in HeLa using TBL1X antibody.

注意事项 (Note)

For research use only .