

产品名称: GTF2H1 (3E8) Mouse Monoclonal Antibody
产品货号: AMM03439



产品概述 (Summary)

产品名称 (Production Name)	GTF2H1 (3E8) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTF2H1
别名 (Alternative Names)	P62; BTF2; TFB1; TFIIH.
基因 ID (Gene ID)	2965
蛋白 ID (SwissProt ID)	P32780.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 62 kDa

研究背景 (Background)

Component of the core-TFIIH basal transcription factor involved in nucleotide excision repair (NER) of DNA and, when

产品名称: GTF2H1 (3E8) Mouse Monoclonal Antibody
产品货号: AMM03439




complexed to CAK, in RNA transcription by RNA polymerase II.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of GTF2H1 in HeLa, HeLa NE and Jurkat lysates using GTF2H1 antibody.



Immunocytochemistry analysis of GTF2H1 (3E8) in HeLa using GTF2H1 antibody.



Immunoprecipitation analysis of GTF2H1 (3E8) in HeLa lysates using GTF2H1 antibody.

产品名称: GTF2H1 (3E8) Mouse Monoclonal Antibody
产品货号: AMM03439



注意事项 (Note)

For research use only .