

产品名称: PPP5C (2F2F7) Mouse Monoclonal Antibody
产品货号: AMM03444



产品概述 (Summary)

产品名称 (Production Name)	PPP5C (2F2F7) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP5C FLJ36922; FLJ55954; OTTHUMP00000165834; PP-T; PP5; PPP5; PPP5_HUMAN;
别名 (Alternative Names)	PPP5C; PPT; Protein phosphatase 5; catalytic subunit; Protein phosphatase T; Serine/threonine protein phosphatase 5; Serine/threonine-protein phosphatase 5.
基因 ID (Gene ID)	5536
蛋白 ID (SwissProt ID)	P53041.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 57 kDa

产品名称: PPP5C (2F2F7) Mouse Monoclonal Antibody
产品货号: AMM03444




研究背景 (Background)

May play a role in the regulation of RNA biogenesis and/or mitosis. In vitro, dephosphorylates serine residues of skeletal muscle phosphorylase and histone H1.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of PP5 in HeLa, Jurkat, 293T, A431 and MCF-7 lysates using PP5 antibody.

注意事项 (Note)

For research use only .