

产品名称: Aurora A (2E3) Mouse Monoclonal Antibody
产品货号: AMM03449



产品概述 (Summary)

产品名称 (Production Name)	Aurora A (2E3) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AURKA
别名 (Alternative Names)	AIK; ARK1; AYK1; Aurora-A; Aurora-related kinase 1; BTAK; IAK1; Ipl1- and aurora-related kinase 1; STK15; STK6; Serine/threonine kinase 15
基因 ID (Gene ID)	6790
蛋白 ID (SwissProt ID)	O14965.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa; Observed MW: 46 kDa

研究背景 (Background)

产品名称: Aurora A (2E3) Mouse Monoclonal Antibody
产品货号: AMM03449




Mitotic serine/threonine kinases that contributes to the regulation of cell cycle progression. Associates with the centrosome and the spindle microtubules during mitosis and plays a critical role in various mitotic events including the establishment of mitotic spindle, centrosome duplication, centrosome separation as well as maturation, chromosomal alignment, spindle assembly checkpoint, and cytokinesis.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of Aurora Kinase A in SW480 and COS7 lysates using Aurora Kinase A antibody.



Immunocytochemistry analysis of Aurora A (2E3) in HeLa using Aurora Kinase A antibody.



产品名称: Aurora A (2E3) Mouse Monoclonal Antibody
产品货号: AMM03449



Immunofluorescence analysis of Aurora A (2E3) in Hela cells using Aurora Kinase A antibody(green) and Histone H3.1 (Phospho-Ser1) antibody (red),and DAPI (blue)

注意事项 (Note)

For research use only .