

产品名称: Syntenin 1 (3D9) Mouse Monoclonal Antibody
产品货号: AMM03460



产品概述 (Summary)

产品名称 (Production Name)	Syntenin 1 (3D9) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SDCBP SDCBP; MDA9; SYCL; Syntenin-1; Melanoma differentiation-associated protein 9; MDA-9; Pro-TGF-alpha cytoplasmic domain-interacting protein 18; TACIP18; Scaffold protein Pbp1; Syndecan-binding protein 1
别名 (Alternative Names)	
基因 ID (Gene ID)	6386
蛋白 ID (SwissProt ID)	O00560.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 32 kDa; Observed MW: 32 kDa

产品名称: Syntenin 1 (3D9) Mouse Monoclonal Antibody
产品货号: AMM03460




研究背景 (Background)

Seems to function as an adapter protein. In adherens junctions may function to couple syndecans to cytoskeletal proteins or signaling components. Seems to couple transcription factor SOX4 to the IL-5 receptor (IL5RA).

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Syntenin in HeLa lysates using Syntenin antibody.

注意事项 (Note)

For research use only .