

产品名称: hnRNP U (1B11) Mouse Monoclonal Antibody
产品货号: AMM03461



产品概述 (Summary)

产品名称 (Production Name)	hnRNP U (1B11) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HNRNPU Heterogeneous nuclear ribonucleoprotein U; hnRNP U protein; HNRNPU;
别名 (Alternative Names)	hnRNPU protein; HNRPU; HNRPU_HUMAN; p120; p120 nuclear protein; SAF A; SAF-A; SAFA; Scaffold attachment factor A; U21.1.
基因 ID (Gene ID)	3192
蛋白 ID (SwissProt ID)	Q00839.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 91 kDa; Observed MW: 110 kDa

产品名称: hnRNP U (1B11) Mouse Monoclonal Antibody
产品货号: AMM03461




研究背景 (Background)

Component of the CRD-mediated complex that promotes MYC mRNA stabilization. Binds to pre-mRNA. Has high affinity for scaffold-attached region (SAR) DNA. Binds to double- and single-stranded DNA and RNA.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of hnRNP U in K562, Jurkat and HeLa lysates using hnRNP U antibody.

注意事项 (Note)

For research use only .