

产品名称: ALDH3A2 (3F7) Mouse Monoclonal Antibody
产品货号: AMM03465



产品概述 (Summary)

产品名称 (Production Name)	ALDH3A2 (3F7) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALDH3A2
别名 (Alternative Names)	ALDH3A2; ALDH10; FALDH; Fatty aldehyde dehydrogenase; Aldehyde dehydrogenase 10; Aldehyde dehydrogenase family 3 member A2; Microsomal aldehyde dehydrogenase
基因 ID (Gene ID)	224
蛋白 ID (SwissProt ID)	P51648.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 55 kDa

产品名称: ALDH3A2 (3F7) Mouse Monoclonal Antibody
产品货号: AMM03465




研究背景 (Background)

Catalyzes the oxidation of long-chain aliphatic aldehydes to fatty acids. Active on a variety of saturated and unsaturated aliphatic aldehydes between 6 and 24 carbons in length.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Aldehyde dehydrogenase 10 in A549, U251 and HepG2 lysates using Aldehyde dehydrogenase 10 antibody.

注意事项 (Note)

For research use only .