

产品名称: Ezrin (2C12) Mouse Monoclonal Antibody
产品货号: AMM03466



产品概述 (Summary)

产品名称 (Production Name)	Ezrin (2C12) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EZR
别名 (Alternative Names)	EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
基因 ID (Gene ID)	7430
蛋白 ID (SwissProt ID)	P15311.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa; Observed MW: 81 kDa

研究背景 (Background)

Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for

产品名称: Ezrin (2C12) Mouse Monoclonal Antibody
产品货号: AMM03466




the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Ezrin in MCF-7, 3T3, A549 and rat Brain lysates using Ezrin antibody.

注意事项 (Note)

For research use only .