

产品名称: PRMT1 (4E9) Mouse Monoclonal Antibody
产品货号: AMM03469



产品概述 (Summary)

产品名称 (Production Name)	PRMT1 (4E9) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRMT1
别名 (Alternative Names)	ANM 1; ANM1; ANM1_HUMAN; complete cds; HCP 1; HCP1; Heterogeneous nuclear ribonucleoprotein methyltransferase 1 like 2; Heterogeneous nuclear ribonucleoproteins methyltransferase like 2; Heterogeneous nuclear ribonucleoproteins methyltransferase like2; Histone-arginine N-methyltransferase PRMT1; HMT 2; HMT1 (hnRNP methyltransferase; HMT1 (hnRNP methyltransferase S. cerevisiae) like 2; HMT1 hnRNP methyltransferase; HMT1 hnRNP methyltransferase like 2 (S. cerevisiae); HMT1 hnRNP methyltransferase like 2; HMT1 hnRNP methyltransferase-like 2 (S. cerevisiae); HMT2; HRMT1 L2; HRMT1L 2; HRMT1L2; Human mRNA for suppressor for yeast mutant; Human mRNA for suppressor for yeast mutant complete cds; Interferon receptor 1 bound protein 4; Interferon receptor 1

产品名称: PRMT1 (4E9) Mouse Monoclonal Antibody
产品货号: AMM03469



bound protein4; Interferon receptor 1-bound protein 4; Interferon receptor 1bound protein 4; IR1 B4; IR1B 4; IR1B4; Mrmt 1; Mrmt1; PRMT 1; PRMT1; Protein arginine methyltransferase 1; Protein arginine N methyltransferase 1; Protein arginine N methyltransferase1; Protein arginine N-methyltransferase 1; R1B4; S. cerevisiae like 2.

基因 ID (Gene ID) 3276
蛋白 ID (SwissProt ID) Q99873.

产品应用 (Application)

稀释比 (Dilution Ratio) WB 1:500-1:1000
蛋白分子量 (Molecular Weight) Calculated MW: 42 kDa; Observed MW: 42 kDa


研究背景 (Background)

Arginine methyltransferase that methylates (mono and asymmetric dimethylation) the guanidino nitrogens of arginyl residues present in proteins such as ESR1, histone H2, H3 and H4, PIAS1, HNRNPA1, HNRNPD, NFATC2IP, SUPT5H, TAF15 and EWS. Constitutes the main enzyme that mediates monomethylation and asymmetric dimethylation of histone H4 'Arg-4' (H4R3me1 and H4R3me2a, respectively), a specific tag for epigenetic transcriptional activation.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of PRMT1 in HeLa, A549 and SW480 lysates using PRMT1 antibody.

注意事项 (Note)

For research use only .