

产品名称: FYB (2H9) Mouse Monoclonal Antibody
产品货号: AMM03470



产品概述 (Summary)

产品名称 (Production Name)	FYB (2H9) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FYB1
别名 (Alternative Names)	ADAP; Fyb; FYN binding protein; p120/p130; SLAP130
基因 ID (Gene ID)	2533
蛋白 ID (SwissProt ID)	O15117.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 85 kDa; Observed MW: 120 kDa

研究背景 (Background)

Acts as an adapter protein of the FYN and LCP2 signaling cascades in T cells. Modulates the expression of interleukin-2 (IL-

产品名称: FYB (2H9) Mouse Monoclonal Antibody
产品货号: AMM03470




2). Involved in platelet activation. Prevents the degradation of SKAP1 and SKAP2.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of FYB in 10, 15 and 20ug Jurkat lysates using FYB antibody.

注意事项 (Note)

For research use only .