

产品名称: Beclin 1 (4F4) Mouse Monoclonal Antibody
产品货号: AMM03473



产品概述 (Summary)

产品名称 (Production Name)	Beclin 1 (4F4) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BECN1
别名 (Alternative Names)	BECN1
基因 ID (Gene ID)	8678
蛋白 ID (SwissProt ID)	Q14457.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 55 kDa

研究背景 (Background)

Plays a central role in autophagy. Required for the abscission step in cytokinesis. May play a role in antiviral host defense.

产品名称: Beclin 1 (4F4) Mouse Monoclonal Antibody
产品货号: AMM03473




Protects against infection by a neurovirulent strain of Sindbis virus.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Beclin 1 in HeLa, 3T3, C6, COS7, CHO-K1, rat brain and Jurkat lysates using Beclin 1 antibody.

注意事项 (Note)

For research use only .