

产品名称: Thioredoxin Reductase 2 (3F2) Mouse Monoclonal Antibody
产品货号: AMM03476



产品概述 (Summary)

产品名称 (Production Name)	Thioredoxin Reductase 2 (3F2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TXNRD2
别名 (Alternative Names)	mitochondrial; selenoprotein Z; SELZ; Thioredoxin reductase 2; Thioredoxin reductase 2 mitochondrial; thioredoxin reductase 3; thioredoxin reductase beta; Thioredoxin reductase TR3; TR 3; TR; TR beta; TR-beta; TRXR 2; TRXR2; TRXR2_HUMAN; TXNRD 2; Txnrd2.
基因 ID (Gene ID)	10587
蛋白 ID (SwissProt ID)	Q9NNW7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
----------------------	-----------------

产品名称: Thioredoxin Reductase 2 (3F2) Mouse Monoclonal Antibody
产品货号: AMM03476



蛋白分子量 (Molecular Weight) Calculated MW: 57 kDa; Observed MW: 57 kDa


研究背景 (Background)

Maintains thioredoxin in a reduced state. Implicated in the defenses against oxidative stress. May play a role in redox-regulated cell signaling.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of TXNRD2 in A431, U251, U2OS, Lncap and Hela lysates using TXNRD2 antibody.

注意事项 (Note)

For research use only .