

产品名称: SMARCC1 (2A4) Mouse Monoclonal Antibody
产品货号: AMM03479



产品概述 (Summary)

产品名称 (Production Name)	SMARCC1 (2A4) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMARCC1
别名 (Alternative Names)	BAF155; CRACC1; Rsc8; SMARCC1; SRG3; SWI3
基因 ID (Gene ID)	6599
蛋白 ID (SwissProt ID)	Q92922.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 123 kDa; Observed MW: 155 kDa

研究背景 (Background)

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-

产品名称: SMARCC1 (2A4) Mouse Monoclonal Antibody
产品货号: AMM03479




nucleosome topology). May stimulate the ATPase activity of the catalytic subunit of the complex.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of SMARCC1 in K562, Jurkat and HeLa lysates using SMARCC1 antibody.

注意事项 (Note)

For research use only .