

**产品名称: Switch Associated Protein 70 (2E9) Mouse
Monoclonal Antibody**
产品货号: AMM03480



产品概述 (Summary)

产品名称 (Production Name)	Switch Associated Protein 70 (2E9) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SWAP70
别名 (Alternative Names)	HSPC321; SWAP-70; FLJ39540; KIAA0640
基因 ID (Gene ID)	23075
蛋白 ID (SwissProt ID)	Q9UH65.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa; Observed MW: 69 kDa

研究背景 (Background)

产品名称: Switch Associated Protein 70 (2E9) Mouse Monoclonal Antibody
产品货号: AMM03480




Phosphatidylinositol 3,4,5-trisphosphate-dependent guanine nucleotide exchange factor (GEF) which, independently of RAS, transduces signals from tyrosine kinase receptors to RAC.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of SWAP70 in COS7, 293T, C6, 3T3 and HeLa lysates using SWAP70 antibody.



Immunoprecipitation analysis of Switch Associated Protein 70 (2E9) in HeLa lysates using SWAP70 antibody.

注意事项 (Note)

For research use only .