

产品名称: Cytochrome P450 17A1 (6E1) Mouse Monoclonal Antibody
产品货号: AMM03481



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome P450 17A1 (6E1) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYP17A1
别名 (Alternative Names)	CPT7; CYP17; P450C17; S17AH; CYP17A1
基因 ID (Gene ID)	1586
蛋白 ID (SwissProt ID)	P05093.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 60 kDa

研究背景 (Background)

产品名称: Cytochrome P450 17A1 (6E1) Mouse Monoclonal Antibody
产品货号: AMM03481




Conversion of pregnenolone and progesterone to their 17-alpha-hydroxylated products and subsequently to dehydroepiandrosterone (DHEA) and androstenedione.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Cytochrome P450 17A1 in Hela, HepG2 and A549 lysates using Cytochrome P450 17A1 antibody.

注意事项 (Note)

For research use only .