

产品名称: CHD4 (5F1) Mouse Monoclonal Antibody
产品货号: AMM03491



产品概述 (Summary)

产品名称 (Production Name)	CHD4 (5F1) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHD4
别名 (Alternative Names)	CHD4; Chromodomain-helicase-DNA-binding protein 4; CHD-4; ATP-dependent helicase CHD4; Mi-2 autoantigen 218 kDa protein; Mi2-beta
基因 ID (Gene ID)	1108
蛋白 ID (SwissProt ID)	Q14839.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 218 kDa; Observed MW: 260 kDa

研究背景 (Background)

产品名称: CHD4 (5F1) Mouse Monoclonal Antibody
产品货号: AMM03491




Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CHD4 in K562, HeLaNE, HL-60, MOLT4, Jurkat and 293 lysates using CHD4 antibody.



Immunoprecipitation analysis of CHD4 (5F1) in K562 lysates using CHD4 antibody.

注意事项 (Note)

For research use only .