

产品名称: Noggin (1D2) Mouse Monoclonal Antibody
产品货号: AMM03498



产品概述 (Summary)

产品名称 (Production Name)	Noggin (1D2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG3
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NOG
别名 (Alternative Names)	SYM1; SYNS1.
基因 ID (Gene ID)	9241
蛋白 ID (SwissProt ID)	Q13253.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

研究背景 (Background)

Essential for cartilage morphogenesis and joint formation. Inhibitor of bone morphogenetic proteins (BMP) signaling which

产品名称: **Noggin (1D2) Mouse Monoclonal Antibody**
产品货号: **AMM03498**




is required for growth and patterning of the neural tube and somite.


研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Western blot analysis of Noggin in CHO-K1 lysates overexpressing NogginPDGFR transmembrane domain fused protein using Noggin antibody.



Immunocytochemistry analysis of Noggin (1D2) in HeLa using Noggin antibody.

注意事项 (Note)

For research use only .