

产品名称: Olig2 (7A1) Mouse Monoclonal Antibody
产品货号: AMM03505



产品概述 (Summary)

产品名称 (Production Name)	Olig2 (7A1) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OLIG2
别名 (Alternative Names)	Basic domain helix loop helix protein class B 1; Basic helix loop helix protein class B 1; BHLHB; bHLHB1; bHLHe19; Class B basic helix loop helix protein 1; Class B basic helix-loop-helix protein 1; class E basic helix loop helix protein 19
基因 ID (Gene ID)	10215
蛋白 ID (SwissProt ID)	Q13516.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 32 kDa; Observed MW: 32 kDa

产品名称: Olig2 (7A1) Mouse Monoclonal Antibody
产品货号: AMM03505




研究背景 (Background)

Olig2 is required for oligodendrocyte and motor neuron specification in the spinal cord, as well as for the development of somatic motor neurons in the hindbrain. Cooperates with OLIG1 to establish the pMN domain of the embryonic neural tube. Antagonist of V2 interneuron and of NKX2-2-induced V3 interneuron development.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of OLIG2 in CHO-K1 lysates and CHO-K1 transfected by EGFPOLIG2 fragment lysates using OLIG2 antibody.

注意事项 (Note)

For research use only .