

产品名称: Cytokeratin pan (7H5) Mouse Monoclonal Antibody
产品货号: AMM03522



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin pan (7H5) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT5
别名 (Alternative Names)	58 kDa cytokeratin; CK 5; CK-5; CK5; Cytokeratin-5; Cytokeratin5; DDD; DDD1; EBS 2; EBS2; epidermolysis bullosa simplex 2 Dowling-Meara/Kobner/Weber-Cockayne types; K2C5_HUMAN; K5
基因 ID (Gene ID)	3852
蛋白 ID (SwissProt ID)	P13647.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 62 kDa

产品名称: Cytokeratin pan (7H5) Mouse Monoclonal Antibody
产品货号: AMM03522




研究背景 (Background)

The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Keratin(5060kDa) in C6, 3T3, A549, MCF-7 and HeLa lysates using Keratin(5060kDa) antibody.

注意事项 (Note)

For research use only .