

产品名称: DDX3 (6G8) Mouse Monoclonal Antibody
产品货号: AMM03525



产品概述 (Summary)

产品名称 (Production Name)	DDX3 (6G8) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, IP, ChIP
种属反应性 (Reactivity)	Human, Rat, Mouse, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDX3X
别名 (Alternative Names)	DBX; DDX3; HLP2; DDX14.
基因 ID (Gene ID)	1654
蛋白 ID (SwissProt ID)	O00571.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50, ChIP 1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 73 kDa; Observed MW: 73 kDa

研究背景 (Background)

Multifunctional ATP-dependent RNA helicase.


产品名称: DDX3 (6G8) Mouse Monoclonal Antibody
产品货号: AMM03525




研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of DDX3 in HeLa, 3T3, C6, COS7, K562 and Jurkat lysates using DDX3 antibody.



Immunocytochemistry analysis of DDX3 (6G8) in HeLa using DDX3 antibody.



Immunoprecipitation analysis of DDX3 (6G8) in HeLa lysates using DDX3 antibody.

注意事项 (Note)

产品名称: DDX3 (6G8) Mouse Monoclonal Antibody
产品货号: AMM03525



For research use only .