

产品名称: BAG1 (2C12) Mouse Monoclonal Antibody
产品货号: AMM03531



产品概述 (Summary)

产品名称 (Production Name)	BAG1 (2C12) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAG1
别名 (Alternative Names)	BAG1; HAP; BAG family molecular chaperone regulator 1; BAG-1; Bcl-2-associated athanogene 1
基因 ID (Gene ID)	573
蛋白 ID (SwissProt ID)	Q99933.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 52,46,33 kDa

研究背景 (Background)

产品名称: BAG1 (2C12) Mouse Monoclonal Antibody
产品货号: AMM03531




Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Inhibits the pro-apoptotic function of PPP1R15A, and has anti-apoptotic activity. Markedly increases the anti-cell death function of BCL2 induced by various stimuli.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of BAG1 (2C12) in CHO-K1, transfected with a Human pEGFPC1BAG1 construct (A) or transfected with a Human pEGFPC1 construct (B), using Bag1 antibody.

注意事项 (Note)

For research use only.