

产品名称: 4E BP2 (7A6) Mouse Monoclonal Antibody
产品货号: AMM03542



产品概述 (Summary)

产品名称 (Production Name)	4E BP2 (7A6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EIF4EBP2
别名 (Alternative Names)	4EBP2; PHASII
基因 ID (Gene ID)	1979
蛋白 ID (SwissProt ID)	Q13542.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 13 kDa; Observed MW: 15-20 kDa

研究背景 (Background)

Regulates eIF4E activity by preventing its assembly into the eIF4F complex. Mediates the regulation of protein translation

产品名称: 4E BP2 (7A6) Mouse Monoclonal Antibody
产品货号: AMM03542




by hormones, growth factors and other stimuli that signal through the MAP kinase pathway.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of 4EBP2 in CHO-K1 lysates and CHO-K1 transfected by EGFP4EBP2 lysates using 4EBP2 antibody.

注意事项 (Note)

For research use only .