

产品名称: Lamin B2 (8F6) Mouse Monoclonal Antibody
产品货号: AMM03554



产品概述 (Summary)

产品名称 (Production Name)	Lamin B2 (8F6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LMNB2
别名 (Alternative Names)	LMN2; LAMB2
基因 ID (Gene ID)	84823
蛋白 ID (SwissProt ID)	Q03252.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

Lamins are components of the nuclear lamina, a fibrous layer on the nucleoplasmic side of the inner nuclear membrane,

产品名称: Lamin B2 (8F6) Mouse Monoclonal Antibody
产品货号: AMM03554




which is thought to provide a framework for the nuclear envelope and may also interact with chromatin.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Lamin B2 in HCT116, mouse Embryo and Jurkat lysates using Lamin B2 antibody.



Immunoprecipitation analysis of Lamin B2 (8F6) in HeLa lysates using Lamin B2 antibody.

注意事项 (Note)

For research use only .