

产品名称: Bid (3C5) Mouse Monoclonal Antibody
产品货号: AMM03565



产品概述 (Summary)

产品名称 (Production Name)	Bid (3C5) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Ascitic Fluid

免疫原信息 (Immunogen)

基因名 (Gene Name)	BID
别名 (Alternative Names)	FP497; MGC15319; MGC42355; BID
基因 ID (Gene ID)	637
蛋白 ID (SwissProt ID)	P55957.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 22 kDa

研究背景 (Background)

The major proteolytic product p15 BID allows the release of cytochrome c (By similarity). Isoform 1, isoform 2 and isoform 4 induce ICE-like proteases and apoptosis. Isoform 3 does not induce apoptosis. Counters the protective effect of Bcl-2.


产品名称: Bid (3C5) Mouse Monoclonal Antibody
产品货号: AMM03565




研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Bid in 293T, Jurkat and A431 lysates using Bid antibody.



Immunocytochemistry analysis of Bid (3C5) in HeLa using Bid antibody.

注意事项 (Note)

For research use only .