

产品名称: LPP (6F6) Mouse Monoclonal Antibody
产品货号: AMM03568



产品概述 (Summary)

产品名称 (Production Name)	LPP (6F6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Ascitic Fluid

免疫原信息 (Immunogen)

基因名 (Gene Name)	LPP
别名 (Alternative Names)	LPP; Lipoma-preferred partner; LIM domain-containing preferred translocation partner in lipoma
基因 ID (Gene ID)	4026
蛋白 ID (SwissProt ID)	Q93052.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 66 kDa; Observed MW: 66 kDa

研究背景 (Background)

May play a structural role at sites of cell adhesion in maintaining cell shape and motility. In addition to these structural

产品名称: LPP (6F6) Mouse Monoclonal Antibody
产品货号: AMM03568




functions, it may also be implicated in signaling events and activation of gene transcription.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of LPP in A549, MCF-7, C6 and HeLa lysates using LPP antibody.




Immunofluorescence analysis of LPP (6F6) in COS using LPP antibody(green). Red: Actin filaments have been labeled using DY554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.




Immunofluorescence analysis of LPP (6F6) in HeLa cells using LPP antibody(green). Red: Actin filaments have been labeled using DY554 phalloidin, and DAPI (blue). DRAQ5 fluorescent DNA dye.


产品名称: LPP (6F6) Mouse Monoclonal Antibody
产品货号: AMM03568



Immunohistochemistry analysis of paraffin-embedded Human vessels tissues using LPP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human small intestine using LPP (6F6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunoprecipitation analysis of LPP in HeLa lysates using LPP (6F6) antibody.

注意事项 (Note)

For research use only .