

产品名称: DiMethyl-Histone H3 (Lys27) (5F6) Mouse Monoclonal Antibody
产品货号: AMM03581



产品概述 (Summary)

产品名称 (Production Name)	DiMethyl-Histone H3 (Lys27) (5F6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Methylated
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H3C1
别名 (Alternative Names)	H3K27me2; H3 histone; HIST1H3A; Histone cluster 1; H3a
基因 ID (Gene ID)	8350
蛋白 ID (SwissProt ID)	P68431.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

研究背景 (Background)

产品名称: DiMethyl-Histone H3 (Lys27) (5F6) Mouse Monoclonal Antibody
产品货号: AMM03581




H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Immunoprecipitation analysis of DiMethyl-Histone H3 (Lys27) (5F6) in HeLa lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody.




Western blot analysis of DiMethyl-Histone H3 (Lys27) (5F6) in C6, 3T3, C2C12 lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody

**产品名称: DiMethyl-Histone H3 (Lys27) (5F6) Mouse
Monoclonal Antibody**
产品货号: AMM03581



Immunocytochemistry analysis of DiMethyl-Histone H3 (Lys27) in HeLa cells using DiMethyl-Histone H3 (Lys27) (5F6) antibody.



Western blot analysis of DiMethyl-Histone H3 (Lys27) (5F6) in HeLa, rat heart, mouse heart, HUVEC, Jurkat lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody.

注意事项 (Note)

For research use only .