

产品名称: Survivin (8B9) Mouse Monoclonal Antibody
产品货号: AMM03600



产品概述 (Summary)

产品名称 (Production Name)	Survivin (8B9) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BIRC5
别名 (Alternative Names)	BIRC5; API4; IAP4; Baculoviral IAP repeat-containing protein 5; Apoptosis inhibitor 4; Apoptosis inhibitor survivin
基因 ID (Gene ID)	332
蛋白 ID (SwissProt ID)	O15392.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 16 kDa

研究背景 (Background)

产品名称: Survivin (8B9) Mouse Monoclonal Antibody
产品货号: AMM03600




Survivin is an apoptosis inhibitor that is expressed during the G2/M phase of the cell cycle. Associates with the microtubules of the mitotic spindle and any disruption results in the loss of apoptosis activity. May play a role in neoplasia.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using Survivin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Survivin (8B9) in TF1, HL6 lysates using Survivin (8B9) antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using Survivin antibody. High-pressure and

产品名称: Survivin (8B9) Mouse Monoclonal Antibody
产品货号: AMM03600



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .