

产品名称: GRP78 BiP (6H7) Mouse Monoclonal Antibody
产品货号: AMM03620



产品概述 (Summary)

产品名称 (Production Name)	GRP78 BiP (6H7) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPA5 HSPA5; GRP78; 78 kDa glucose-regulated protein; GRP-78; Endoplasmic
别名 (Alternative Names)	reticulum lumenal Ca(2+)-binding protein grp78; Heat shock 70 kDa protein 5; Immunoglobulin heavy chain-binding protein; BiP
基因 ID (Gene ID)	3309
蛋白 ID (SwissProt ID)	P11021.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 72 kDa; Observed MW: 78 kDa

产品名称: GRP78 BiP (6H7) Mouse Monoclonal Antibody
产品货号: AMM03620




研究背景 (Background)

When Chinese hamster K12 cells are starved of glucose, the synthesis of several proteins, called glucose-regulated proteins (GRPs), is markedly increased. Hendershot et al. (1994) (PubMed 8020977) pointed out that one of these, GRP78 (HSPA5), also referred to as 'immunoglobulin heavy chain-binding protein' (BiP), is a member of the heat-shock protein-70 (HSP70) family and is involved in the folding and assembly of proteins in the endoplasmic reticulum (ER).


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma using GRP78 BiP (6H7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemical analysis of paraffin-embedded Human tonsils using GRP78 BiP (6H7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.


产品名称: GRP78 BiP (6H7) Mouse Monoclonal Antibody
产品货号: AMM03620




Western blot analysis of GRP78 BiP (6H7) in Jurkat, MCF-7, SKBR3, HeLa lysates using GRP78 BiP (6H7) antibody



Immunocytochemistry analysis of GRP78 BiP in HeLa cells using GRP78 BiP (6H7) antibody.



Immunoprecipitation analysis of GRP78 BiP in HeLa lysates using BiGRP78 BiP (6H7) antibody.



产品名称: GRP78 BiP (6H7) Mouse Monoclonal Antibody
产品货号: AMM03620



Western blot analysis of GRP78 BiP (6H7) in CHO-K1, rat brain, C6, 3T3, Hela lysates using GRP78 BiP (6H7) antibody

注意事项 (Note)

For research use only .